


INSTRUCTIONS FOR USE



 **aokinmycontrol** **OTA IAC**

Order No.: MY-IC-M-04

Sample preparation with aokinImmunoClean**M** columns
(IAC)

Analytical-kit for rapid and quantitative determination of Ochratoxin (OTA).

Materials

aokinmycontrolOTA IAC (Order No.: MY-IC-M-04-100)

Package content

A) *IC Consumables:*

Filter paper

aokinExtractionSalt OTA + spoon

aokinICelute OTA

Reaction tubes 2 mL

IC-Reservoirs (loading reservoirs) 25 mL

B) *IC Liquids:*

aokinExtractionSolvent OTA

aokinICDilute 10x

aokinICWash

C) *Immunoaffinity columns:*

aokinImmunoClean M OTA (Order No.: IC-M-04-100)

D) *Materials for analytical measurement:*

aokinReactionBuffer, Reaction buffer

aokinmycontrolOTA/AFLA SPE, Additive (black cap),

ADD-OTA/AFLA SPE, (for 5 analyses each)

aokinmycontrolOTA, Reagent 1 (yellow cap), F-OTA,

(for 5 analyses each)

aokinmycontrolOTA, Reagent 2 (white cap), A-OTA,

(for 5 analyses each)




Figure 1: Reagent 1, F-OTA (yellow cap), Reagent 2, A-OTA (white cap) and Reaction buffer (1 L bottle)

E) *Materials for internal quality control:*

aokinmycontrol OTA, negative control OTA

(transparent), (for zero value measurements)

aokinmycontrolOTA/AFLA SPE, Additive (black cap),

ADD-OTA/AFLA SPE, (for 5 analyses each)

aokinmycontrolOTA, Reagent 1 (yellow cap), F-OTA,

(for 5 analyses each)

aokinmycontrolOTA, Reagent 2 (white cap), A-OTA,

(for 5 analyses each)

Note: All substances provided are precisely weighed and calibrated. Control of the volume and concentration of the individual solutions are essential for the precision of the analysis.

Caution: The extraction solvent may contain methanol. Work with professional care.

Storage Conditions: Reagents 1, 2 and aokinImmunoClean M OTA columns must be stored at temperature of +2 – +10°C. All other components may be stored at room temperature.

Quality Control: All materials and reagents are prepared according to strict quality control protocols. Exchanging reagents between kits having different Lot-numbers will lead to erroneous results and is not permitted.

Order Information:

aokinmycontrolOTA (Order No.: MY-IC-M-04-100)

Introduction

aokinmycontrolOTA is a rapid and precise quantitative method for analyzing Ochratoxin (OTA). It has been specifically designed and calibrated for use with aokinImmunoClean M columns. The sample preparation can be completed in 25-45 minutes.

aokinmycontrolOTA is available with a calibration. Please use professional care and check the accuracy by regularly analyzing reference materials (aokinReference-MatrixMaterials) and/or standards. Participation in proficiency tests is recommended.

aokin will gladly assist you to customise this test for your specific sample type and application. Please do not hesitate to contact us.

Sample	corn, grain, coffee
Time required for sample preparation	25 to 45 minutes
Time required for measurement	3 minutes
Analysis	
	Measurement range [µg/kg]
Range 1	1-10
Range 2	2-20
Range 3	4-40

Ochratoxin

Ochratoxin is a mycotoxin produced by fungal molds such as *Aspergillus* and *Penicillium*. Ochratoxin is known to be commonly present in commodities such as cereals, coffee, dried fruit and red wine. It can accumulate in the body and is considered a human carcinogen. Exposure to ochratoxins through diet can have acute toxicity to kidneys. It is strongly recommended to monitor ochratoxin contamination in food and feed products.




Figure 2: Chemical formula for Ochratoxin A, C₂₀H₁₈ClNO₆. Molecular weight: 403.81g/mol

Recommended Accessories

All required materials are available from aokin.

Tel.: +49 30 9489 2160

	Order No.
aokinextractor (food blender)	EX-07-06
aokinwatchbox (timer for food blender)	EX-07-06-4
Weighing scale, d = 0,01 g	LB-03-04
Eppendorf centrifuge, variable g-force	LB-04-04
Variable pipettes (1000 µl)	LB-04-05-1000
Pipette tips (1000 µl)	LB-04-08-1000
Funnel	LB-05-04
Dispensette	LB-08-01
aokinReferenceMatrixMaterial	RMM-04
Reaction tubes, 2mL	LB-05-02
aokinICadapter	LB-08-10-05
aokinVacuumAdapter	LB-08-11
Vacuum-pump (diaphragm pump)	LB-04-10
Trap for vacuum-pump (vacuum-bottle)	LB-04-12
Vacuum manifold	LB-04-09
Glass fiber filters, GF/F	LB-04-13-100

Vacuum bottle and funnel for glass fiber filters (Büchner-Funnel)	LB-04-14
pH-indicator strips	LB-08-18
Magnetic stirrer plate (not heated)	LB-04-15
Magnetic stirrer (to degas)	LB-04-16-05
aokinCWashTween	TW-05-1000
PBS	PBS-01-1000

Sample preparation

The following protocol is an example. The quantification ranges are dependent on dilutions. Actual volume settings in the software may vary.

Note: It is of critical importance to use the correct sample preparation protocol for each quantification range. Use volumes displayed in the *aokin* software.

1. Quality control

There are free materials included in the kit, for your internal quality control: **Reagent 1**, **Reagent 2**, **Additive**, as well as **negative control** solutions for measurements of zero values (corresponding to samples free of mycotoxin).

Please perform measurements of negative controls regularly, this ensures the accuracy of your determinations.

If you notice increased values, change cuvette and repeat measurement. If sample results remain high, contact the *aokin* team.

2. Sample collection, homogenisation, and grinding

The analysis sample is collected, ground, and homo-genised according to an approved procedure. Small samples may be ground using the *aokinextractor*.

3. Weighing and extraction

Weigh 15 g of your sample, add one spoon (1,5 g) of *aokinExtractionSalt* OTA and 20 g extraction solution (22 ml *aokinExtractionSolvent* OTA at 20°C) respectively 40g (44mL) for absorbent samples like cocoa directly into the extraction beaker.



Figure 3: Weighing

Close the extraction beaker with the lid (with the blending knives). Blend for 3.5 minutes. The recommended protocol has blending times alternating with resting time to avoid heating of the sample and is as follows: mix for 30 seconds, pause for 1 minute, mix for 30 seconds and so on (until 3.5 minutes of blending time).

Use the *aokinwatchbox* (a preprogrammed timer) to conveniently and automatically complete this extraction protocol.



Figure 4: Extracting with the aokinextractor (blender)

4. Dilution

Add 166 mL *aokinCDilute10x* directly into extraction beaker. Shake softly.

5. Filtration

Place the filter on a suitable funnel and the funnel onto a collection container. Open the extraction beaker, pour the extract onto the filter and collect the filtrate (at least 25 mL). Discard the filter paper and filter cake. Shake/stir the filtrate to ensure homogeneity.



Figure 5: Filtration

Check the pH value. If necessary, readjust the pH value by adding diluted HCl or NaOH to a neutral pH (pH 6.5 - 7.5).

Important: all precipitate has to be removed before loading the column by glass fiber filtration to ensure a good flow through the column (see step 5b).



Figure 6: Filtration through a glass fiber filter using a vacuum (shown as an example for a coffee sample)

5. Clean up with *aokinImmunoClean M*

a) Preparation

Attach the vacuum manifold to a pump. Attach the *aokinImmunoClean M* column to the vacuum manifold (or alternatively onto an *aokin VacuumAdapter*) and transfer *aokinCWash* solution onto the column. Attach the *aokinCadapter* to the *aokinImmunoClean M* column and

connect it to the loading reservoir. Pipette aokinCWash into the loading reservoir.

Important: No air bubbles should be visible in the column. Let the aokinCWash solution elute, leaving only a little liquid in the reservoir and so that the aokinCadapter and aokinImmunoClean M column are still filled with liquid.

b) Loading

Load the column with 20 ml of the diluted filtrate. Adjust the vacuum so that a flow rate of approximately **1 ml/min** is produced (total running time approximately 20 minutes). The running speed can be reduced by closing the stopcock, or increased with applying vacuum. Alternatively, a positive pressure from above can be used.

Important: Do not let the aokinImmunoClean M column run dry before, or during loading.

Important: High flow rates lead to a reduced rate of recovery.

c) Washing

Fill 2 x 5 ml aokinCWashTween into the loading reservoir. Adjust the vacuum so that a high flow rate is achieved (running time approximately 2 minutes). The column bed should appear relatively colorless. If this is not observed, the wash step may be repeated. Complete the wash step by adding 2 x 5mL PBS twice to the column. Respectively wash with 2 x 15mL IC-wash, see page 6 till 9. Let the column run dry.



Figure 7: Attachment of reservoir, aokinCadapter and aokinImmunoClean M column to the aokinVacuumAdapter (left), to the vacuum manifold with a coffee sample (middle), and the coffee column after washing with aokinCWash or aokinCWashTween and PBS (right).

d) Elution

Place the column into a collection tube and centrifuge for 1 minute at 1000 x g to remove the residual liquid.



Figure 8: The column is placed into a collection tube before eluent is added, and then the column is centrifuged to remove residual liquid.

Set the column into a new clean collection tube and pipette 300 µl aokinCElute OTA into the column, close (half closed) the column after that with the lid. When closing the column, about 100 µl aokinCElute OTA will be **pushed through the gel**. After incubation for 3 minutes, centrifuge for 1 minute at 3000 x g. The lid must remain on the column during centrifugation. The liquid is present in the collection vessel after centrifugation.

Important note: Eluates of green coffee samples and others may show a precipitate. This leads to increased results. If you see or suspect a precipitation in the eluate than centrifuge eluate at high speed, than pipette supernatant into new collection tube and use for further analysis.

6. Analyzing

Use column-filtrate for analyzing in the aokinspectrometerFP470.

Important: Gas bubbles can lead to measurement inaccuracies. It is strongly recommended to degas the reaction buffer (aokinReactionBuffer) before each use. This is performed by applying vacuum and stirring for 30 min.



Figure 9: Setup for degassing the reaction buffer (aokinReactionBuffer). This consists of: vacuum pump (hose and vacuum trap, not shown) and the reaction buffer in a bottle with a magnetic stirrer rod on a magnetic stirrer plate.

Please follow detailed instructions for spectrometer use (aokinspectrometerFP470 & aokinLHW03 Instructions for use).

This includes:

- 1) Place *Reagents 1* and *2* into position A6 and B6, and place **Additive** into position B1 of the sample rack of your spectrometer.
- 2) Place a 25 ml *Clean1*-glass container, filled with *Clean1* solution into the *Clean1*-position, on the left side, next to the palette.
- 3) Place a 25 ml *neg. control OTA*-glass container, filled with *negative control OTA* solution = aokinCElute OTA in the *Clean2*-position, on the left side of the palette.
- 4) Place an empty 2 mL vial in position A1 of the palette.
- 5) Place an empty waste bottle into the holder. Check presence of Reaction buffer and check if tubing is below the surface.

-
- 6) Place a new cuvette with a clean stirrer into the spectrometer.

7. Quality control

Included in the analytical kit there are following additional materials for your internal quality control: Reagent 1, Reagent 2, negative control samples (labelled negative control, corresponding to samples free of mycotoxin) and a positive control sample.

Please perform measurements of negative controls regularly, this ensures the accuracy of your determinations.

If you notice increased values, change cuvette and repeat measurement. If results remain high perform an offset correction of the calibration based on the negative control results. In addition, the use of recovery corrections preferentially by using sample extracts from certified reference matrix samples or alternatively based on the positive controls included in the kit.

If problem persists calibrate. Please contact the aokin team for any support needed.

Conversion factor: analyte concentration in cuvette (nM) to amount in sample (µg/kg)

 **aokinmycontrol** **OTA IAC**
Standard

Step 1: Extraction

- Sample mass: $m_{\text{Sample}} = 15 \text{ g}$
- Volume extraction solvent: $V_{\text{Extraction solvent}} = 188 \text{ ml}$
(22 ml_{EL OTA} + 16 ml_{CDilute 10x} + 150 ml_{Water})
- Molar mass Ochratoxin: $MW_{\text{OTA}} = 403,81 \left[\frac{\text{g}}{\text{mol}} \right]$

Mycotoxin concentration in the sample extract:

$$c \left[\frac{\mu\text{mol}}{\text{l}} \right]_{\text{Extract}} = \frac{m_{\text{Sample}} [\text{kg}]}{V_{\text{Solvent}} [\text{l}] * MW_{\text{Mycotoxin}} \left[\frac{\text{g}}{\text{mol}} \right]} * c \left[\frac{\mu\text{g}}{\text{kg}} \right]_{\text{Sample}} = \frac{0.015}{0.188 * 403.81} * c \left[\frac{\mu\text{g}}{\text{kg}} \right]_{\text{Sample}} = 0.00019759 * c \left[\frac{\mu\text{g}}{\text{kg}} \right]_{\text{Sample}}$$

Step 2: Purification
with **aokin IC M OTA**

- Volume sample extract load to the **aokin IC-M** column: $V_{\text{loaded sample extract}} = 20 \text{ ml}$
- Volume "aokin Elute" used for the elution of the column: $V_{\text{elute}} = 0,3 \text{ ml}$

Mycotoxin concentration in the column eluate:

$$c \left[\frac{\mu\text{mol}}{\text{l}} \right]_{\text{Eluate}} = \frac{V_{\text{load}} [\text{ml}]}{V_{\text{elute}} [\text{ml}]} * c \left[\frac{\mu\text{mol}}{\text{l}} \right]_{\text{Extract}} = \frac{20}{0.3} * c \left[\frac{\mu\text{mol}}{\text{l}} \right]_{\text{Extract}} = 66,66 * c \left[\frac{\mu\text{mol}}{\text{l}} \right]_{\text{Extract}}$$

Step 3: Measurement
aokin FP 470 / LHW 03

- Sample volume: $V_{\text{Column eluate}} = V_{\text{Sample}} = 200 \mu\text{l}$
- Total volume in the cuvette: $V_{\text{Cuvette}} = 2600 \mu\text{l}$

Mycotoxin concentration in the cuvette:

$$c \left[\frac{\mu\text{mol}}{\text{l}} \right]_{\text{Cuvette}} = \frac{V_{\text{Sample}} [\mu\text{l}]}{V_{\text{Cuvette}} [\mu\text{l}]} * c \left[\frac{\mu\text{mol}}{\text{l}} \right]_{\text{Eluate}} = \frac{200}{2600} * c \left[\frac{\mu\text{mol}}{\text{l}} \right]_{\text{Eluate}} = 0.0769 * c \left[\frac{\mu\text{mol}}{\text{l}} \right]_{\text{Eluate}}$$


Conversion factor: Extraction, Purification and Measurement
It follows the conversion factor from 1 to 3 above:

$$c \left[\frac{\mu\text{mol}}{\text{l}} \right]_{\text{Cuvette}} = 0.00019759 * 66,66 * 0.0769 * c \left[\frac{\mu\text{g}}{\text{kg}} \right]_{\text{Sample}} = 0.001013 * c \left[\frac{\mu\text{g}}{\text{kg}} \right]_{\text{Sample}} \quad \text{or}$$







$$c \left[\frac{\mu\text{g}}{\text{kg}} \right]_{\text{Sample}} = \frac{1}{1.013} * c \left[\frac{\text{nmol}}{\text{l}} \right]_{\text{Cuvette}} = 0.987 * c \left[\frac{\text{nmol}}{\text{l}} \right]_{\text{Cuvette}}$$







Ochratoxin / standard samples:

- Recommended for standard types of samples:
wheat, corn, feeds, nuts, paprika, green and roasted coffee
- Not suitable for absorbent samples → please use the sample preparation for **absorbent samples (aokinmycontrol OTA Absorbent)**
- Not suitable for colored samples → please use the sample preparation for **special matrices (aokinmycontrol OTA Special Matrix)**

 **aokinmycontrol** **OTA**
Standard

Procedure:


Extraction		Weighing: 15 g sample 1,5 g aokinExtractionSalt OTA 22 mL aokinExtractionSolvent OTA
		Extraction: 3,5 min mixing with aokinwatchbox
		166 ml aokinICDilute 10x
		Filtration: collect filtrate, (discard filter cake)
		check and adjust to 6.5 - 7.5 pH, neutralize if necessary by adding NaOH or HCl
		filtrate through a glass fiber filter (optional)

Purification		set up column
		20 ml filtrate, 1 drop / second (slow flow rate)
		2 x 15 ml aokinICWash
		1 min centrifuge at 1000 x g
		<p>Elution:</p> <p>300 µl aokinICElute OTA</p> <p>3 min incubation with lid on</p> <p>1 min centrifuge at 3000 x g, use eluate for measurement</p>
Measurement		<p>Automatic Analyse (FP470 / LHW03):</p> <p>place the 2ml reaction tube in the sample holder of the <i>LHW03</i></p> <p>2000 µl aokin Reaction buffer</p> <p>200 µl aokinmycontrol OTA Additive</p> <p>200 µl sample</p> <p style="text-align: right;"><i>(diluted 1:1 - RANGE 1)</i> <i>(diluted 1:2 - RANGE 2)</i> <i>(diluted 1:4 - RANGE 3)</i></p> <p>100 µl aokinmycontrol OTA Reagent 1</p> <p>100 µl aokinmycontrol OTA Reagent 2</p>







OTA = Ochratoxin, Conversion factor: 1 nmol OTA/l in cuvette = 0,987 µg/kg in sample







Ochratoxin / absorbent samples:

- Recommended for absorbent samples like cocoa
- Not suitable for standard samples → please use the sample preparation for **standard samples (aokinmycontrol OTA Standard)**
- Not suitable for colored samples → please use the sample preparation for **special matrices (aokinmycontrol OTA Special Matrix)**

 **aokinmycontrol** **OTA**
Absorbent

Procedure:

Extraction		Weighing: 15 g sample 1,5 g aokinExtractionSalt OTA 44 mL aokinExtractionSolvent OTA
		Extraction: 3,5 min mixing with aokinwatchbox
		166 ml aokinICDilute 10x
		Filtration: collect filtrate, (discard filter cake)
		check and adjust to 6.5 - 7.5 pH, neutralize if necessary by adding NaOH or HCl
		filtrate through a glass fiber filter (optional)

Purification		set up column
		20 ml filtrate, 1 drop / second (slow flow rate)
		2 x 15 ml aokinCWash
		1 min centrifuge at 1000 x g
		Elution: 300 µl aokinCElute OTA 3 min incubation with lid on 1 min centrifuge at 3000 x g, use eluate for measurement
Measurement		Automatic Analyse (FP470 / LHW03): place the 2ml reaction tube in the sample holder of the LHW03 2000 µl aokin Reaction buffer 200 µl aokinmycontrol OTA Additive 200 µl sample (diluted 1:1 - RANGE 1) (diluted 1:2 - RANGE 2) (diluted 1:4 - RANGE 3) 100 µl aokinmycontrol OTA Reagent 1 100 µl aokinmycontrol OTA Reagent 2







OTA = Ochratoxin, Conversion factor: 1 nmol OTA/l in cuvette = 1,974 µg/kg in sample







Ochratoxin / special matrix:

- Recommended for colored sample types like spices (except paprika), some feeds and some nuts
- Not suitable for standard samples → please use the sample preparation for **standard samples (aokinmycontrol OTA Standard)**
- Not suitable for absorbent samples → please use the sample preparation for **absorbent samples (aokinmycontrol OTA Absorbent)**

 **aokinmycontrol** **OTA**
Special Matrix

Procedure:

Extraction		Weighing: 15 g sample 1,5 g aokinExtractionSalt OTA 22 mL aokinExtractionSolvent OTA
		Extraction: 3,5 min mixing with aokinwatchbox
		166 ml aokinICDilute 10x
		Filtration: collect filtrate, at least 25 ml (discard filter cake)
		check and adjust to 6.5 - 7.5 pH, neutralize if necessary by adding NaOH or HCl
		filtrate through a glass fiber filter (optional)

Purification		set up column
		20 ml filtrate, 1 drop / second (slow flow rate)
		5 - 20 ml aokinICWashTween 2 x 5 ml PBS
		1 min centrifuge at 1000 x g
		Elution: 300 µl aokinICElute OTA 3 min incubation with lid on 1 min centrifuge at 3000 x g, use eluate for measurement
Measurement		Automatic Analyse (FP470 / LHW03): place the 2ml reaction tube in the sample holder of the <i>LHW03</i> 2000 µl aokin Reaction buffer 200 µl aokinmycontrol OTA Additive 200 µl sample (diluted 1:1 - RANGE 1) (diluted 1:2 - RANGE 2) (diluted 1:4 - RANGE 3) 100 µl aokinmycontrol OTA Reagent 1 100 µl aokinmycontrol OTA Reagent 2

OTA = Ochratoxin, Conversion factor: 1 nmol OTA/l in cuvette = 0,987 µg/kg in sample