


ImmunoAffinity Column for Ochratoxin A for sample preparation

Order no.: IC-C+-04-25

Lot: xxx xxx xxx xxx



aokinImmunoClean C+ **OTA** column



Ochratoxin A (OTA)

Specification:

Recommended loading: < 200 ng OTA
Loading rate: < 3 ml/min
Elution rate: < 1 ml/min

	Recovery of 150 ng
Specification:	<input checked="" type="checkbox"/> Recovery > 80%

Methodology Testing:

Sample loading: OTA standard 0 ng, 1 ng, 150 ng Ochratoxin A
Fluid 20% Methanol / 10mM PBS pH 7,2 / 0,05% Sodiumazide
Volume 16 ml
Rate 1 ml / min

Column washing: Fluid 10mM PBS, pH 7,2 / 0,05% Sodiumazide
Volume 20 ml
Rate > 1 ml /min

Elution: Eluant 2% Glacial acetic acid/ 98% Methanol
Volume 1.5 ml
Rate 1 ml / min
Incubation 2 min

Analysis: by FPIA



Dipl.-Ing. Biotechnology
Nicole Abdallah
Head of Quality Control