


## ImmunoAffinity Column for Ochratoxin A for sample preparation

**Order no.: IC-C-04-25**  
Lot: xxx xxx xxx xxx



**aokinImmunoClean C **OTA** column**



Ochratoxin A (OTA)

### **Specification:**

Recommended loading: < 200 ng OTA  
Loading rate: ~ 1 ml/min  
Elution rate: < 1 ml/min

	Recovery of 150 ng
Specification:	<input checked="" type="checkbox"/> Recovery > 80%

### **Methodology Testing:**

**Sample loading:** OTA standard 0 ng, 150 ng Ochratoxin A (OTA)  
Fluid 10 % Methanol, Phosphate Buffer, pH 7.4  
Volume 10 ml  
Rate 1 ml / min

**Column washing:** Fluid 10 % Methanol, Phosphate Buffer, pH 7.4  
Volume 5 ml  
Rate > 1 ml / min

**Elution:** Eluant 0,004% Glacial Acetic Acid / Methanol  
Volume 2 ml  
Incubation 2 min  
Rate 1 ml / min

**Analysis:** by FPIA



Dipl.-Ing. Biotechnology  
Nicole Abdallah  
Head of Quality Control