




## ImmunoAffinity Column for T2/HT2 toxin for sample preparation


**Order no.: IC-C-78-25**  
Lot: xxx xxx xxx xxx



aokinImmunoClean C T2/HT2 column



T2-Toxin (T2)



HT2-Toxin (HT2)

### Specification:

Recommended loading: < 500 ng T2/HT2  
Loading rate: 1-2 ml/min  
Elution rate: < 1 ml/min

	Recovery of 100 ng	Recovery of 500 ng
Specification:	<input checked="" type="checkbox"/> Recovery > 80%	<input checked="" type="checkbox"/> Recovery > 80%

### Methodology Testing:

**Sample loading:**

T2 standard	100 ng, 500 ng T2 standard
Fluid	10% Methanol / 10 mM PBS, pH 7.2 / 0,05% Sodiumazide
Volume	10 ml
Rate	1 ml / min

**Column washing:**

Fluid	10 mM PBS
Volume	5 ml
Rate	> 1 ml / min

**Elution:**

Eluant	100% Methanol
Volume	1.5 ml
Incubation	2 min
Rate	1 ml / min

**Analysis:** by HPLC



Dipl.-Ing. Biotechnology  
Nicole Abdallah  
Head of Quality Control