

Order No.: CJ-05-BSA Lot: 005XXXXX305

## **Conjugate FUM-BSA**

Conjugate of Fumonisin B1 (FUM) with Bovine Serum Albumin (BSA)

## Specification

Protein	Bovine Serum Albumin (BSA)
Label	Fumonisin B1
Molar Excess of NHS ester during coupling	50-fold
Protein Concentration	$cm \approx 1.0 mg/ml$
Volume	V = 500 µl
Amount of Conjugate	$m \approx 0.5 \text{ mg}$
Buffer	PBS pH 7.4 without preservatives
Quality Controller	Lot specific data
Storage	4 °C, protect from light
Retest date	Lot specific data

## Note

If the final concentration of purified protein conjugate is less than 1 mg/ml, we recommend to add bovine serum albumin (BSA) or other stabilizing proteins at 1-10 mg/ml. In the presence of 0.02 % thimerosal or 2 mM sodium azide, the conjugate should be stable at 4 °C for several months. For longer storage, divide the conjugate into small aliquots and freeze at -20 °C. **Avoid repeated freezing and thawing!** 

## Caution

The use of this product is strictly limited to trained personnel for professional manufacturing, laboratory, or research purposes. Substance not yet fully tested. Final Fitness-For-Use must be determined by and is the sole responsibility of the end-user.

aokin AG  $\cdot$  Robert-Rössle-Str. 10  $\cdot$  13125 Berlin  $\cdot$  Germany

Tel. +49 30 94 89 21 60 · Fax +49 30 94 89 21 61 · info@aokin.com · www.aokin.com