

## Conjugate T2-BSA

Conjugate of T2-Toxin (T2) with Bovine Serum Albumin (BSA)

### Specification

Protein	Bovine Serum Albumin (BSA)
Label	T2-Toxin
Molar Excess of NHS ester during coupling	20-fold
Protein Concentration	cm $\approx$ 1.0 mg/ml
Volume	V = 500 $\mu$ l
Amount of Conjugate	m $\approx$ 0.5 mg
Buffer	PBS pH 7.4 without preservatives
Quality Controller	Lot specific data
Storage	4 °C, protect from light
Retest date	Lot specific data

### Note

If the final concentration of purified protein conjugate is less than 1 mg/ml, we recommend to add bovine serum albumin (BSA) or other stabilizing proteins at 1-10 mg/ml. In the presence of 0.02 % thimerosal or 2 mM sodium azide, the conjugate should be stable at 4 °C for several months. For longer storage, divide the conjugate into small aliquots and freeze at -20 °C. **Avoid repeated freezing and thawing!**

### Caution

The use of this product is strictly limited to trained personnel for professional manufacturing, laboratory, or research purposes. Substance not yet fully tested. Final Fitness-For-Use must be determined by and is the sole responsibility of the end-user.