

aokin ImmunoClean (IC) columns are designed for specific isolation and concentration of the analyte prior to analysis. The purified sample may then be analyzed using *aokin* mycontrol FP, HPLC, or LC-MS, or another analytical method. *aokin* ImmunoClean columns are available as either gravity flow or spin columns.

aokin QuickClean (QC) columns rapidly improve the quality of sample extracts for analysis. They eliminate matrix effects which can cause false negatives or false positives and ensure the reproducibility of any analytical method. *aokin* QuickClean columns are easy to use. Just add sample and centrifuge for 2 minutes immediately prior to detection with an *aokin* analyzerFP, HPLC, LC-MS or other detection system. Total hands-on time is under 3 minutes.

Product code:

... C: gravity flow column

... M: spin column, a small centrifuge capable of 5000 G is required

Products

Mycotoxin	Type	Size	Units/pack	Order No.
Aflatoxin B1, B2, G1, G2	ImmunoClean	C	25	IC-C-03-25
		M	50	IC-M-03-50
	QuickClean	M	100	QC-03-100
Aflatoxin M1	ImmunoClean	C	25	IC-C-31-25
		M	50	IC-M-31-50
Aflatoxin total / Ochratoxin	ImmunoClean	C	25	IC-C-0304-25
		M	50	IC-M-0304-50
Deoxynivalenol	ImmunoClean	C	25	IC-C-02-25
		M	50	IC-M-02-50
	QuickClean	M	100	QC-02-100
Fumonisin	ImmunoClean	C	25	IC-C-05-25
		M	50	IC-M-05-50
Ochratoxin	ImmunoClean	C	25	IC-C-04-25
		M	50	IC-M-04-50
Ochratoxin / Aflatoxin total	ImmunoClean	C	25	IC-C-0304-25
		M	50	IC-M-0304-50
Paraquat / Diquat	QuickClean	M	100	QC-70-100
Patulin	QuickClean	M	100	QC-06-100
T-2 / HT-2 Toxin	ImmunoClean	C	25	IC-C-78-25
		M	50	IC-M-78-50
	QuickClean	M	100	QC-78-100

Mycotoxin	Type	Size	Units/pack	Order No.
Zearalenone	ImmunoClean	C	25	IC-C-01-25
		M	50	IC-M-01-50
	QuickClean	M	100	QC-01-100
Zearanol	ImmunoClean	C	25	IC-C-11-25
		M	50	IC-M-11-50
