

aokin mycontrol analyzer *FP*

Rapid Accurate Detection

A portable, hand-held, high-sensitivity instrument for the rapid, specific, and quantitative analysis of mycotoxins, herbicides, THC and other analytes.

Simple, field-based analysis with HPLC-grade accuracy in just 3 minutes.

- Low cost per sample
- Built-in standard curves and calibration
- High performance – excellent correlation with HPLC



*aokin*mycontrol analyzer fluorescence polarometer

aokin mycontrol detection kits

Order No.	Name	Assay for the detection of	Number of tests
MYS-03	aokin mycontrol AFLA	Aflatoxin (total)	20/100/500
MYS-031	aokin mycontrol AFLA M1	Aflatoxin M1	20/100/500
MYS-02	aokin mycontrol DON	Deoxynivalenol	20/100/500
MYS-05	aokin mycontrol FUM	Fumonisin	20/100/500
MYS-99	aokin mycontrol IGG Quant	Immunoglobulin G	20/100/500
MYS-78	aokin mycontrol T2/HT2	T2/HT2 Toxin	20/100/500
MYS-40	aokin mycontrol THC	Tetrahydrocannabinol (total)	20/100/500
MYS-01	aokin mycontrol ZON	Zearalenone	20/100/500
MYS-11	aokin mycontrol ZEL	Zeranol (α -zearalanol, zearalanol)	20/100/500
MYS-15	aokin mycontrol CAF	Caffeine	20/100/500
	aokin mycontrol CBD	in development	



aokin AG · Robert-Rössle-Str. 10 · 13125 Berlin · Germany
 Tel. +49 30 94 89 21 60 · Fax +49 30 94 89 21 61 · info@aokin.com