




Conjugate FUM-BSA

Order-No: CJ-05-BSA

Lot: xxx xxx xxx xxx

Structure:



Conjugate of Fumonisin B1 (FUM) with Bovine Serum Albumin (BSA)

Specification:

- Synthesis number: xxx xxx xxx xxx
- Protein: Bovine Serum Albumin (BSA)
- Label: Fumonisin B1
- Molar Excess of NHS ester during coupling: 50-fold
- Protein Concentration: $c_m \approx 1.0$ mg/ml
- Volume: $V = 200$ μ l
- Amount of Conjugate: $m \approx 0.2$ mg
- Buffer: PBS pH 7.4 without preservatives
- Quality Controller: D.C.
- Storage: 4° C, protect from light!
- Expiry date: 1 year after delivery
- Hint: If the final concentration of purified protein conjugate is less than 1 mg/ml, we recommend to add bovine serum albumin (BSA) or other stabilizing proteins at 1^{-10} mg/ml. In the presence of 0.02% thimerosal or 2 mM sodium azide, the conjugate should be stable at 4 °C for several months. For longer storage, divide the conjugate into small aliquots and freeze at -20 °C. **Avoid repeated freezing and thawing!**

Caution

Substance not yet fully tested (EU).
Not for drug, household or other uses.

Aokin AG - 13125 Berlin
Tel: +49 (0) 3094892160