


## T2/HT2 Toxin (T2/HT2)

HT2-Toxin: C<sub>22</sub>H<sub>32</sub>O<sub>8</sub>, Molecular weight: 434, 48 g/mol

<b>aokin Standard solutions</b>	
<b>aokin</b> Standard <b>T2</b> , 1 ml (100 µg/ml) in Acetonitrile	Order-No.: <b>CH-07-L1-100</b>
<b>aokin</b> Standard <b>T2</b> , 5 ml (100 µg/ml) in Acetonitrile	Order-No.: <b>CH-07-L5-100</b>
<b>aokin</b> Standard <b>HT2</b> , 1 ml (100 µg/ml) in Acetonitrile	Order-No.: <b>CH-08-L1-100</b>
<b>aokin</b> Standard <b>HT2</b> , 5 ml (100 µg/ml) in Acetonitrile	Order-No.: <b>CH-08-L5-100</b>

<b>aokin Standards (solid)</b>	
<b>aokin</b> Standard <b>T2</b> , 1 mg	Order-No.: <b>CH-07-S1</b>
<b>aokin</b> Standard <b>T2</b> , 5 mg	Order-No.: <b>CH-07-S5</b>
<b>aokin</b> Standard <b>HT2</b> , 1 mg	Order-No.: <b>CH-08-S1</b>
<b>aokin</b> Standard <b>HT2</b> , 5 mg	Order-No.: <b>CH-08-S5</b>

<b>aokin Antibodies</b>	
<b>aokin</b> Monoclonal Antibody <b>T2</b> (Clone2), 200 µg, in PBS	Order-No.: <b>AB-07-02</b>
<b>aokin</b> Monoclonal Antibody <b>HT2</b> (Clone1), 100% cross reactive to T2, 200 µl, in TRIS-Citrate	Order-No.: <b>AB-08-01</b>




## T2/HT2 Toxin (T2/HT2)

HT2-Toxin: C<sub>22</sub>H<sub>32</sub>O<sub>8</sub>, Molecular weight: 434, 48 g/mol

<b>aokin</b> Conjugates	
<b>aokin</b> Conjugate <b>T2-BSA (Bovine Serum Albumin)</b> , 500 µl (1 mg/ml) in PBS	Order-No.: <b>CJ-07-BSA</b>
<b>aokin</b> Conjugate <b>T2-HRP (Horseradish Peroxidase)</b> , 500 µl (1 mg/ml) in PBS	Order-No.: <b>CJ-07-HRP</b>
<b>aokin</b> Conjugate <b>T2-KLH (Keyhole Limpet Hemocyanine)</b> , 500 µl (1 mg/ml) in PBS	Order-No.: <b>CJ-07-KLH</b>
<b>aokin</b> Conjugate <b>T2-OVA (Ovalbumin)</b> , 500 µl (1 mg/ml) in PBS	Order-No.: <b>CJ-07-OVA</b>
<b>aokin</b> Conjugate <b>T2-RPE (R-Phycoerythrin)</b> , 500 µl (1 mg/ml) in PBS	Order-No.: <b>CJ-07-RPE</b>
<b>aokin</b> Conjugate <b>HT2-BSA (Bovine Serum Albumin)</b> , 500 µl (1 mg/ml) in PBS	Order-No.: <b>CJ-08-BSA</b>
<b>aokin</b> Biotinylated <b>T2/HT2 Toxins</b> and Intermediates ( <b>Oximes</b> )	on demand

<b>aokin</b> ImmunoClean columns	
<b>aokin</b> ImmunoClean C <b>T2/HT2</b> (25 units / pack)	Order-No.: <b>IC-C-78-25</b>
<b>aokin</b> ImmunoClean C+ <b>T2/HT2</b> (25 units / pack)	Order-No.: <b>IC-C+-78-25</b>
<b>aokin</b> ImmunoClean CF <b>T2/HT2</b> (25 units / pack)	Order-No.: <b>IC-CF-78-25</b>
<b>aokin</b> ImmunoClean M <b>T2/HT2</b> (50 units / pack)	Order-No.: <b>IC-M-78-50</b>



## T2/HT2 Toxin (T2/HT2)

HT2-Toxin: C<sub>22</sub>H<sub>32</sub>O<sub>8</sub>, Molecular weight: 434, 48 g/mol

### **aokin** QuickClean columns

**aokin** QuickClean **T2/HT2** (100 units / pack)

Order-No.: **QC-78-100**

### **aokin** mycontrol

**aokin** mycontrol **T2/HT2** (20 measurements)

Order-No.: **MY-QC-78-20**

**aokin** mycontrol **T2/HT2** (100 measurements)

Order-No.: **MY-QC-78-100**

### **aokin** ReferenceMatrixMaterials

[\*\*aokin\*\* ReferenceMatrixMaterials \*\*T2/HT2\*\*](#)

Order-No.: **RMM-78**