

## Zearalenone (ZON)


Zearalenone: C<sub>18</sub>H<sub>22</sub>O<sub>5</sub>, Molecular weight: 318,36 g/mol

<b>aokin Standard solutions</b>	
<b>aokin</b> Standard <b>ZON</b> , 1 ml (100 µg/ml) in Acetonitrile	Order-No.: <b>CH-01-L1-100</b>
<b>aokin</b> Standard <b>ZON</b> , 5 ml (100 µg/ml) in Acetonitrile	Order-No.: <b>CH-01-L5-100</b>

<b>aokin Standards (solid)</b>	
<b>aokin</b> Standard <b>ZON</b> , 1 mg	Order-No.: <b>CH-01-S1</b>
<b>aokin</b> Standard <b>ZON</b> , 5 mg	Order-No.: <b>CH-01-S5</b>

<b>aokin Antibodies</b>	
<b>aokin</b> Monoclonal Antibody <b>ZON</b> (Clone1), 200 µg in PBS	Order-No.: <b>AB-01-01</b>

<b>aokin Conjugates</b>	
<b>aokin</b> Conjugate <b>ZON-BSA (Bovine Serum Albumin)</b> , 500 µl (1 mg/ml) in PBS	Order-No.: <b>CJ-01-BSA</b>
<b>aokin</b> Conjugate <b>ZON-HRP (Horseradish Peroxidase)</b> , 500 µl (1 mg/ml) in PBS	Order-No.: <b>CJ-01-HRP</b>
<b>aokin</b> Conjugate <b>ZON B1-KLH (Keyhole Limpet Hemocyanine)</b> , 500 µl (1 mg/ml) in PBS	Order-No.: <b>CJ-01-KLH</b>
<b>aokin</b> Conjugate <b>ZON-OVA (Ovalbumin)</b> , 500 µl (1 mg/ml) in PBS	Order-No.: <b>CJ-01-OVA</b>
<b>aokin</b> Conjugate <b>ZON-RPE (R-Phycoerythrin)</b> , 500 µl (1 mg/ml) in PBS	Order-No.: <b>CJ-01-RPE</b>



## Zearalenone (ZON)

Zearalenone: C<sub>18</sub>H<sub>22</sub>O<sub>5</sub>, Molecular weight: 318,36 g/mol

<b>aokin</b> Biotinylated Zearalenone and Intermediates (Oximes)	on demand
--	-----------

<b>aokin</b> ImmunoClean columns	
<b>aokin</b> ImmunoClean C ZON (25 units / pack)	Order-No.: <b>IC-C-01-25</b>
<b>aokin</b> ImmunoClean C+ ZON (25 units / pack)	Order-No.: <b>IC-C+-01-25</b>
<b>aokin</b> ImmunoClean CF ZON (25 units / pack)	Order-No.: <b>IC-CF-01-25</b>
<b>aokin</b> ImmunoClean M ZON (50 units / pack)	Order-No.: <b>IC-M-01-50</b>

<b>aokin</b> QuickClean columns	
<b>aokin</b> QuickClean ZON (100 units / pack)	Order-No.: <b>QC-01-100</b>

<b>aokin</b> mycontrol	
<b>aokin</b> mycontrol ZON (20 measurements)	Order-No.: <b>MY-QC-01-20</b>
<b>aokin</b> mycontrol ZON (100 measurements)	Order-No.: <b>MY-QC-01-100</b>

<b>aokin</b> ReferenceMatrixMaterials	
<b>aokin</b> ReferenceMatrixMaterials ZON	Order-No.: <b>RMM-01</b>