




**Conjugate T2-BSA**

T2-toxin, conjugate of T2-toxin with Bovine Serum Albumin (BSA)

**Order-No: CJ-07-BSA**

**Lot: 007XXXXXXXXXX305**

**Structure:**



Conjugate of T2-toxin with Bovine Serum Albumin

**Specification:**

- Synthesis number: XXX XX XXX
- Protein: Bovine Serum Albumin (BSA)
- Label: HS-T2
- Bovine Serum Albumin: 20-fold
- Concentration of protein:  $c_m \approx 1.0$  mg/ml
- Amount of Conjugate:  $m \approx 0.5$  mg
- Volume:  $V = 500$   $\mu$ l
- Buffer: PBS pH 7.4, without preservatives
- Storage: 4° C, protect from light!
- Expiry date: 1 year after delivery

**Hint:** Store the labelled protein at 4 °C, protected from light! If the final concentration of purified protein conjugate is less than 1 mg/ml, add bovine serum albumin (BSA) or other stabilizing proteins at 1-10 mg/ml. In the presence of 0.02% thimerosal or 2 mM sodium azide, the conjugate should be stable at 4 °C for several months. For longer storage, divide the conjugate into small aliquots and freeze at - 20 °C. **Avoid repeated freezing and thawing!**

**Caution**

Substance not yet fully tested (EU).  
Not for drug, household or other uses.

Aokin AG – 13125 Berlin  
Tel: +49 (0) 3094892160