




Conjugate HT2-BSA

T2-toxin, conjugate of T2-toxin with Bovine Serum Albumin (BSA)

Order-No: CJ-08-BSA

Lot: 008XXXXXXXXXX305

Structure:



Conjugate of HT2-toxin with Bovine Serum Albumin

Specification:

- Synthesis number: XXX XX XXX
- Protein: Bovine Serum Albumin (BSA)
- Label: HS-T2
- Bovine Serum Albumin: 50-fold
- Concentration of protein: $c_m \approx 1.0$ mg/ml
- Amount of Conjugate: $m \approx 0.5$ mg
- Volume: $V = 500$ μ l
- Buffer: PBS pH 7.4, without preservatives
- Storage: 4° C, protect from light!
- Quality Controller: R.W.
- Expiry date: 1 year after delivery

Hint: Store the labelled protein at 4 °C, protected from light! If the final concentration of purified protein conjugate is less than 1 mg/ml, add bovine serum albumin (BSA) or other stabilizing proteins at 1-10 mg/ml. In the presence of 0.02% thimerosal or 2 mM sodium azide, the conjugate should be stable at 4 °C for several months. For longer storage, divide the conjugate into small aliquots and freeze at - 20 °C. **Avoid repeated freezing and thawing!**

Caution

Substance not yet fully tested (EU).
Not for drug, household or other uses.

Aokin AG - 13125 Berlin
Tel: +49 (0) 3094892160