


T2/HT2 Toxin (T2/HT2)

HT2-Toxin: C₂₂H₃₂O₈, Molecular weight: 434, 48 g/mol

aokin Standard solutions	
aokin Standard T2, 1 ml (100 µg/ml) in Acetonitrile	Order-No.: CH-07-L1-100
aokin Standard T2, 5 ml (100 µg/ml) in Acetonitrile	Order-No.: CH-07-L5-100
aokin Standard HT2, 1 ml (100 µg/ml) in Acetonitrile	Order-No.: CH-08-L1-100
aokin Standard HT2, 5 ml (100 µg/ml) in Acetonitrile	Order-No.: CH-08-L5-100

aokin Standards (solid)	
aokin Standard T2, 1 mg	Order-No.: CH-07-S1
aokin Standard T2, 5 mg	Order-No.: CH-07-S5
aokin Standard HT2, 1 mg	Order-No.: CH-08-S1
aokin Standard HT2, 5 mg	Order-No.: CH-08-S5


aokin Antibodies	
aokin Monoclonal Antibody T2 (Clone2), Anti-T2-Toxin, 200 µg	Order-No.: AB-07-02-200
aokin Monoclonal Antibody T2 (Clone2), Anti-T2-Toxin, 1 mg	Order-No.: AB-07-02-1000
aokin Monoclonal Antibody HT2 (Clone1), Anti-HT2-Toxin, 100% cross reactive to T2, 200 µg	Order-No.: AB-08-01-200
aokin Monoclonal Antibody HT2 (Clone1), Anti-HT2-Toxin, 100% cross reactive to T2, 1 mg	Order-No.: AB-08-01-1000



T2/HT2 Toxin (T2/HT2)

HT2-Toxin: C₂₂H₃₂O₈, Molecular weight: 434, 48 g/mol

aokin Conjugates	
<u>aokin Conjugate T2-BSA (Bovine Serum Albumin), 500 µl (~1 mg/ml)</u>	Order-No.: CJ-07-BSA
<u>aokin Conjugate T2-HRP (Horseradish Peroxidase), 500 µl (~1 mg/ml)</u>	Order-No.: CJ-07-HRP
<u>aokin Conjugate T2-KLH (Keyhole Limpet Hemocyanine), 500 µl (~1 mg/ml)</u>	Order-No.: CJ-07-KLH
<u>aokin Conjugate T2-OVA (Ovalbumin), 500 µl (~1 mg/ml)</u>	Order-No.: CJ-07-OVA
<u>aokin Conjugate T2-RPE (R-Phycoerythrin), 500 µl (~1 mg/ml)</u>	Order-No.: CJ-07-RPE
<u>aokin Conjugate HT2-BSA (Bovine Serum Albumin), 500 µl (~1 mg/ml)</u>	Order-No.: CJ-08-BSA
<u>aokin Conjugate HT2-HRP (Horseradish Peroxidase), 500 µl (~1 mg/ml)</u>	Order-No.: CJ-08-HRP
<u>aokin Conjugate HT2-KLH (Keyhole Limpet Hemocyanine), 500 µl (~1 mg/ml)</u>	Order-No.: CJ-08-KLH
<u>aokin Conjugate HT2-OVA (Ovalbumin), 500 µl (~1 mg/ml)</u>	Order-No.: CJ-08-OVA
<u>aokin Conjugate HT2-RPE (R-Phycoerythrin), 500 µl (~1 mg/ml)</u>	Order-No.: CJ-08-RPE
aokin Biotinylated T2/HT2 Toxins and Intermediates (Oximes)	on demand



T2/HT2 Toxin (T2/HT2)

HT2-Toxin: C₂₂H₃₂O₈, Molecular weight: 434, 48 g/mol

aokin ImmunoClean columns	
aokin ImmunoClean C T2/HT2 (25 units / pack)	Order-No.: IC-C-78-25
aokin ImmunoClean M T2/HT2 (50 units / pack)	Order-No.: IC-M-78-50

aokin QuickClean columns	
aokin QuickClean T2/HT2 (100 units / pack)	Order-No.: QC-78-100

aokin mycontrol	
aokin mycontrol T2/HT2 (20 measurements)	Order-No.: MYS-QC-78-20
aokin mycontrol T2/HT2 (100 measurements)	Order-No.: MYS-QC-78-100

aokin ReferenceMatrixMaterials	
aokin ReferenceMatrixMaterials T2/HT2	Order-No.: RMM-78