

Zearalenone (ZON)


Zearalenone: C₁₈H₂₂O₅, Molecular weight: 318,36 g/mol

aokin Standard solutions	
aokin Standard ZON, 1 ml (100 µg/ml) in Acetonitrile	Order-No.: CH-01-L1-100
aokin Standard ZON, 5 ml (100 µg/ml) in Acetonitrile	Order-No.: CH-01-L5-100

aokin Standards (solid)	
aokin Standard ZON, 1 mg	Order-No.: CH-01-S1
aokin Standard ZON, 5 mg	Order-No.: CH-01-S5

aokin Antibodies	
aokin Monoclonal Antibody ZON (Clone1), Anti-Zearalenone, 200 µg	Order-No.: AB-01-01-200
aokin Monoclonal Antibody ZON (Clone1), Anti-Zearalenone, 1 mg	Order-No.: AB-01-01-1000

aokin Conjugates	
aokin Conjugate ZON-BSA (Bovine Serum Albumin), 500 µl (~1 mg/ml)	Order-No.: CJ-01-BSA
aokin Conjugate ZON-rCTB (recombinant Cholera Toxin B), 500 µl (~1 mg/ml)	Order-No.: CJ-01-CTB
aokin Conjugate ZON-HRP (Horseradish Peroxidase), 500 µl (~1 mg/ml)	Order-No.: CJ-01-HRP
aokin Conjugate ZON-KLH (Keyhole Limpet Hemocyanine), 500 µl (~1 mg/ml)	Order-No.: CJ-01-KLH



Zearalenone (ZON)

Zearalenone: C₁₈H₂₂O₅, Molecular weight: 318,36 g/mol

aokin Conjugate ZON-OVA (Ovalbumin), 500 µl (~1 mg/ml) in PBS	Order-No.: CJ-01-OVA
aokin Conjugate ZON-RPE (R-Phycoerythrin), 500 µl (~1 mg/ml)	Order-No.: CJ-01-RPE
aokin Biotinylated Zearalenone and Intermediates (Oximes)	on demand

aokin ImmunoClean columns	
aokin ImmunoClean C ZON (25 units / pack)	Order-No.: IC-C-01-25
aokin ImmunoClean M ZON (50 units / pack)	Order-No.: IC-M-01-50

aokin QuickClean columns	
aokin QuickClean ZON (100 units / pack)	Order-No.: QC-01-100

aokin mycontrol	
aokin mycontrol ZON (20 measurements)	Order-No.: MYS-QC-01-20
aokin mycontrol ZON (100 measurements)	Order-No.: MYS-QC-01-100

aokin ReferenceMatrixMaterials	
aokin ReferenceMatrixMaterials ZON	Order-No.: RMM-01