


Deoxynivalenol (DON)

Deoxynivalenol: C₁₅H₂₀O₆, Molecular weight: 296,32 g/mol

aokin Standard solutions	
aokin Standard DON , 1 ml (100 µg/ml) in Acetonitrile	Order-No.: CH-02-L1-100
aokin Standard DON , 5 ml (100 µg/ml) in Acetonitrile	Order-No.: CH-02-L5-100

aokin Standards (solid)	
aokin Standard DON , 1 mg	Order-No.: CH-02-S1
aokin Standard DON , 5 mg	Order-No.: CH-02-S5

aokin Antibodies	
aokin Monoclonal Antibody DON (Clone2), 200 µg in PBS	Order-No.: AB-02-02
aokin Monoclonal Antibody DON (Clone4), 200 µg, Protein A purified in TRIS-Citrate	Order-No.: AB-02-04




Deoxynivalenol (DON)

Deoxynivalenol: C₁₅H₂₀O₆, Molecular weight: 296,32 g/mol

aokin Conjugates	
aokin Conjugate DON-BSA (Bovine Serum Albumin) , 500 µl (1 mg/ml) in PBS	Order-No.: CJ-02-BSA
aokin Conjugate DON-HRP (Horseradish Peroxidase) , 500 µl (1 mg/ml) in PBS	Order-No.: CJ-02-HRP
aokin Conjugate DON B1-KLH (Keyhole Limpet Hemocyanine) , 500 µl (1 mg/ml) in PBS	Order-No.: CJ-02-KLH
aokin Conjugate DON-OVA (Ovalbumin) , 500 µl (1 mg/ml) in PBS	Order-No.: CJ-02-OVA
aokin Conjugate DON-TG (Thyroglobulin) , 500 µl (1 mg/ml) in PBS	Order-No.: CJ-02-TG
aokin Conjugate DON-F (Fluorescein) , 500 µl	Order-No.: CJ-02-F
aokin Biotinylated Deoxanivalenol and Intermediates (Oximes)	on demand

aokin ImmunoClean columns	
aokin ImmunoClean C DON (25 units / pack)	Order-No.: IC-C-02-25
aokin ImmunoClean C+ DON (25 units / pack)	Order-No.: IC-C+-02-25
aokin ImmunoClean CF DON (25 units / pack)	Order-No.: IC-CF-02-25
aokin ImmunoClean M DON (50 units / pack)	Order-No.: IC-M-02-50



Deoxynivalenol (DON)

Deoxynivalenol: C₁₅H₂₀O₆, Molecular weight: 296,32 g/mol

aokin QuickClean columns	
aokin QuickClean DON (100 units / pack)	Order-No.: QC-02-100
aokin mycontrol	
aokin mycontrol DON (20 measurements)	Order-No.: MY-QC-02-20
aokin mycontrol DON (100 measurements)	Order-No.: MY-QC-02-100
aokin ReferenceMatrixMaterials	
aokin ReferenceMatrixMaterials DON	Order-No.: RMM-02